# Softwarehandbuch

*Software manual*

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A large blue graphic element consisting of a horizontal bar at the bottom, a diagonal line rising from the left, a horizontal bar at the top, and a diagonal line falling from the top to the right.

# **Software manual**

## **Flexible gripper EGH for cobot**

### **SCHUNK software module for TMflow**

Translation of original software manual

## Imprint

### Copyright:

This manual is protected by copyright. The author is SCHUNK SE & Co. KG.  
All rights reserved.

### Technical changes:

We reserve the right to make alterations for the purpose of technical improvement.

**Document number:** 1449062

**Version:** 07.00 | 14/08/2023 | en

Dear Customer,

Thank you for trusting our products and our family-owned company, the leading technology supplier of robots and production machines.

Our team is always available to answer any questions on this product and other solutions. Ask us questions and challenge us. We will find a solution!

Best regards,

Your SCHUNK team

Customer Management

Tel. +49-7133-103-2503

Fax +49-7133-103-2189

cmg@de.schunk.com

**Please read the operating manual in full and keep it close to the product.**

## Table of Contents

<b>1 General .....</b>	<b>39</b>
1.1 About this manual.....	39
1.2 Presentation of Warning Labels .....	39
1.3 Symbol definition.....	40
1.4 Applicable documents .....	40
<b>2 Description of function .....</b>	<b>41</b>
<b>3 Installing components .....</b>	<b>42</b>
3.1 Preparing the USB stick.....	42
3.2 Importing components .....	43
3.3 Activating components .....	45
3.4 Deactivating components.....	46
<b>4 Uninstalling components.....</b>	<b>47</b>
4.1 Removing a component.....	47
4.2 Removing IODD files .....	48
<b>5 Configuring components .....</b>	<b>49</b>
5.1 Configuring the IP address .....	49
5.2 Connecting LED light bands .....	50
5.3 Setting the Tool Center Point (TCP).....	52
5.4 Optional: Setting functions to "Gripper Button" .....	54
<b>6 Using components in the program .....</b>	<b>55</b>
6.1 Creating the robot program .....	55
6.2 Description of EGH components .....	58
6.2.1 EGH STOP, FSTOP, REF, ACK, STROK.....	58
6.2.2 EGH gripper – V100_GRIP .....	59
6.2.3 EGH Release – V100_RELEASE .....	60
6.2.4 EGH positioning movement – V100_POS .....	60
6.2.5 EGH light band – V100_LIGHT.....	61
6.3 Using multiple grippers .....	62
<b>7 Example programs.....</b>	<b>65</b>
7.1 Example: Initializing gripper.....	65
7.2 Example: Using multiple grippers .....	66

# 1 General

## 1.1 About this manual

This manual contains information on EGH software modules, hereinafter referred to as "components", and their use in "TMflow".

The components are used to easily integrate and control the following products in Techman applications:

- Flexible gripper EGH for cobot

This manual describes the software environment for a Techman robot.

**NOTE:** The illustrations in this manual are intended to provide a basic understanding and may deviate from the actual version.

In addition to these instructions, the documents listed under ▶ 1.4 [ 40] are applicable.

## 1.2 Presentation of Warning Labels

To make risks clear, the following signal words and symbols are used for safety notes.

### **⚠ DANGER**

**Dangers for persons!**

Non-observance will inevitably cause irreversible injury or death.

### **⚠ WARNING**

**Dangers for persons!**

Non-observance can lead to irreversible injury and even death.

### **⚠ CAUTION**

**Dangers for persons!**

Non-observance can cause minor injuries.

### **CAUTION**

**Material damage!**

Information about avoiding material damage.

### 1.3 Symbol definition

The following symbols are used in this manual:

■ Prerequisite for an action

1. Action 1

2. Action 2

⇒ Intermediate results

⇒ Final results

▶ 1.3 [📄 40]: chapter number and [page number] in hyperlinks

### 1.4 Applicable documents

- Assembly and operating manual for the product \*
- Software guide "SCHUNK gripper with IO-Link" \*
- Operating manual for Techman robots

The documents labeled with an asterisk (\*) can be downloaded from [schunk.com/egh-downloads](https://www.schunk.com/egh-downloads).

## 2 Description of function

The components are used for the simple commissioning and programming of the SCHUNK gripper EGH in combination with robots from Techman. All components fit seamlessly into the Techman "TMflow" programming environment. The programming and configuration of the gripper are supported.

### 3 Installing components

Before installing the software module, check the robot software and update it if necessary. The robot software must be compatible with the software of the module.

For information on mutually compatible software versions, see [schunk.com/egh-downloads](http://schunk.com/egh-downloads).

#### CAUTION

##### Possible damage to product!

The product or the robot may get damaged if electrical cables are connected or disconnected during operation.

- Connect or disconnect electrical connections only when the device is switched off.


### 3.1 Preparing the USB stick

#### NOTE

To install the software, SCHUNK recommends using a USB stick.

The USB stick must meet the following requirements:

- Formatted in FAT32 format
- Designation of the removable drive: "TMROBOT".

*Designation USB stick*

#### Storage location for components

1. Download the latest version of the software module at [schunk.com/egh-downloads](http://schunk.com/egh-downloads).
2. Save components in a directory "TM\_Export\Schunk\ComponentObject" on the USB stick.


### Storage location for IODD file (optional)

Note: The IODD file is only required if several grippers are to be operated simultaneously.

1. Download and unpack the latest version of the IODD file at [schunk.com/egh-downloads](http://schunk.com/egh-downloads).
2. Copy the file "SCHUNK-EGH-yyyymmdd-IODD-v.v.xml" into a directory "TM\_Export/Schunk/XmlFiles/IODD" on the USB stick.


## 3.2 Importing components

1. Insert the USB stick into an available USB slot in the control box.
2. Select *Menu – System – Import/Export*.

Call up the "Import/Export" window

3. Select the "Import" button.
  - ⇒ The "Robot List" window opens.
4. Select "SCHUNK" and confirm with "OK".
5. Select the "Component" menu.
6. Select the components to be imported and click the "Import" button.


"Import/Export" window - Importing components

⇒ Components are imported.

### Importing the IODD file (optional)

IODD files only need to be imported if several grippers are going to be operated.

1. Select the menu "IODD Files" in the "Import/Export" window.
2. Select the IODD file to be imported and click the "Import" button.


"Import/Export" window - Import IODD

⇒ IODD file is imported.

### 3.3 Activating components


1. Select *Menu – Parameter – Components*.

Call up the "Component" window

2. Activate the desired component in the "Enable" column.
3. Select the "Save" button.

⇒ The selection box lights up green when the component is activated.


Activating components

⇒ All activated components are displayed in the "Project"  window below.


### 3.4 Deactivating components

1. Select *Menu – Parameter – Components*.

Call up the "Component" window

2. Deactivate the desired component in the "Enable" column.
3. Select the "Save" button.

Deactivating components


⇒ Components are deactivated.

## 4 Uninstalling components

### 4.1 Removing a component


Components can be permanently removed from the program.

1. Select *Menu – Parameter – Components*.

Call up the "Component" window

2. Select the desired component in the right column.
3. Select the "Save" button.
4. Repeat the procedure for other components.


Removing components

⇒ Components are removed permanently.


## 4.2 Removing IODD files

1. Select *Menu – Parameter – Text File Manager*.

Calling up the "Text File Manager" window

2. Select the "IODD Files" tab.
3. Select the desired IODD file and press the red X.

Removing IODD files


⇒ The IODD file is removed.

## 5 Configuring components

### 5.1 Configuring the IP address


For successful communication, the robot must be in the same subnet as the IO-Link master. For this, it is necessary to give the robot a static IP address.

1. Select *Menu – System – Network*.

Call up the "Network Settings" window

2. Connect the IO-Link master to an available Ethernet port on the control box.
3. Select the network adapter to which the device is connected (in the example "Local Area Connection 9").
4. Activate the static IP.
5. Enter the IP address for the robot in the address range **192.168.1.XXX**.
  - ⇒ The default IP address of the IO-Link master is **192.168.1.253**. For both devices to be on the same network, the robot must have an IP address in the address range **192.168.1.XXX**. The last position can be chosen freely.
  - ⇒ Note: The same IP address may not be used by two devices on the network simultaneously. In the example, **192.168.1.156** is used.
6. Enter subnet mask **255.255.255.0**.
7. Confirm entries with "OK".

Network settings

⇒ The robot has been assigned a static IP address.

## 5.2 Connecting LED light bands


If the light band is to be used by the gripper, the cables must be connected to the control box and the outputs from the robot must be switched to PNP.

### Switching ports to PNP

#### NOTE

- In the standard configuration, NPN is the connection (bridge) between GND and 24V. For the switchover to PNP, the bridge (1) between 24V and DOCOM must be connected.
- This change applies to all 16 ports on the control box. That is why only PNP-wired devices work on the control box.

1. **IMPORTANT! Material damage due to short circuit! Do not connect GND to 24V.** Connect the DOCOM connection (1) to the positive terminal of the current supply (24V).
2. The inputs can also be optionally switched to PNP. For this, connect DICOM to GND.

Switching ports to PNP

### Connecting light band cables


By default, the blue cable of the gripper is connected to D00 and the red cable to D01 on the control box.

If other ports are used, they must be changed in the component at runtime, see ▶ [6.2.5 \[ 61 \]](#).

### Connecting the IO-Link master to the control box

- Connect the IO-Link master to the control box via an Ethernet cable. Any available RJ45 socket can be selected.

**IMPORTANT! Do not connect the device to the EtherCat connector of the control box.**


Control box connections

### 5.3 Setting the Tool Center Point (TCP)


For proper use of the gripper on a TM robot, SCHUNK recommends storing the Tool Center Point (TCP) and the gripper weight in the robot settings.

1. Select *Menu – Parameter – TCP Settings*.


Call up the "TCP Settings" window

2. Select the "Manual input parameters of TPC" button.

"TCP Settings" window – Parameters

3. Enter the values in the input fields, see the following table.
4. Select the "Save As" button.
  - ⇒ An input field opens
5. Enter a description for the gripper, e.g. "SCHUNK\_EGH80".
6. Select the "Apply" button to confirm all entries.

TCP and weight input mask

⇒ The Tool Center Point and the gripper weight have been saved in the robot settings.

Parameter	Value during assembly	
	without adapter plate	with adapter plate
X (mm)	0	0
Y (mm)	0	-25
Z (mm)	124.8 + length of gripper fingers	137.3 + length of gripper fingers
RX (deg)	0	0
RY (deg)	0	0
RZ (deg)	0	0
Mass (kg)	0.95	0.95

Tab.: Input values: TCP and weight

## 5.4 Optional: Setting functions to "Gripper Button"

The "Grip" and "Release" functions can be placed on a "Gripper Button". This makes it possible to insert the components in the

"Project"  window by pressing a button on the robot.


Note: The gripper will *not* move mechanically when a button is pressed. This is only possible when the program is running.

1. Select *Menu – Parameter – Gripper Button*.

Call up the "Gripper Button" window

2. Activate "Using Customized Component".
3. Select "Gripper\_Schunk\_V100\_GRIP\_EGH80" for Grip and "Gripper\_Schunk\_V100\_RELEASE\_EGH80" for Release.

Setting functions to "Gripper Button"

⇒ The "Grip" and "Release" functions are placed on a button.

## 6 Using components in the program

### NOTE

In order to transfer the relevant IO-Link data to the IO-Link master, the commands Acknowledge (V100\_ACK) and Reference (V100\_REF) must be set before the first gripping movement.

After this gripping movement, the red "ERROR" LED on the IO-Link master goes out. The gripper is ready for operation.

### 6.1 Creating the robot program

After installation, the components are displayed in TMflow and can be easily dragged and dropped into the program sequence. This creates a sequence plan that the robot processes sequentially.

Additional parameters can be set for the gripping, positioning and light band components. To do so, select the component and choose the "Edit" symbol.

#### 1. Initializing


To start the program, the gripper requires the following clearly defined sequence for initialization:

1. Acknowledge (V100\_ACK)
2. Referencing (V100\_REF)
3. optional: stroke measurement (V100\_STROK)

For an example of how to initialize the robot, see section [7.1 \[ 65 \]](#).









#### 2. Using components

After initialization, the components can be used. These are displayed in the window on the left.

Components are displayed

The following components can be used:


Description	Description of the component in TMflow	Components in TMflow
Stop	V100_STOP	
Quick stop	V100_FSTOP	
Referencing	V100_REF	
Acknowledging	V100_ACK	
Stroke measurement	V100_STROKE	
EGH gripping	V100_GRIP	
EGH release	V100_RELEASE	
EGH positioning movement	V100_POS	
EGH light band	V100_LIGHT	

## Warning for already existing components

### NOTE

A warning appears if a component has been dragged into a sequence that already exists, e.g. a second GRIP component.

1. *When using a single gripper:* Select the "Inherit old component" button.
  - ⇒ The component is overwritten.
2. *When using multiple grippers:* Select the "Add" button.
  - ⇒ For further information, see section ▶ 6.3 [ 62].

"Add new device" warning

## 6.2 Description of EGH components

### 6.2.1 EGH STOP, FSTOP, REF, ACK, STROK

Component	Description
Stop V100_STOP	The gripper is brought to a controlled standstill. The gripping force is maintained.
Quick stop V100_FSTOP	The electrical power supply is interrupted immediately, the gripper is stopped uncontrolled. An error message requiring acknowledgment is generated. <b>IMPORTANT! Material damage to the workpiece possible! Ensure that all workpieces have been removed and that the base jaws are free to move in the reference direction up to the end stop.</b>
Acknowledging V100_ACK	If an error has been eliminated, the gripper is switched from the error state to the normal operating state. The actuator remains de-energized until the next command. In the event of an error that needs to be acknowledged, the first of the three indicator LEDs on the gripper changes from green to red.
Referencing V100_REF	The gripper moves in the configured referencing direction up to the mechanical end stop. The referencing must be run through once at the beginning.  Note: Referencing the product takes some time. During this time, do not perform any other operations on the product or robot.
Stroke measurement V100_STROK	A stroke measurement should be performed if the stroke of the base jaws is limited, for instance, by specific gripper finger shapes. The stroke measurement moves in the opposite direction to the referencing. If no stroke measurement is performed, a stored default value for a stroke of 40 mm per base jaw is used.

## 6.2.2 EGH gripper - V100\_GRIP

With the "Gripping" action, the set gripping direction is moved to the stop and the workpiece is held. Before reaching the entered target position, the motor current and thus the speed of the base jaws is reduced. This protects the gripped workpiece from damage.

### NOTE

During the first gripping process, the complete stroke is driven at a reduced speed and current. The position of the last successful grip is accepted as the intended workpiece position for the next grip command. The intended workpiece position is reset to blockage after referencing, stroke measurement as well as after positioning.

### Parameter settings V100\_GRIP

Gripping parameters	Value	Description
SetDirection	TRUE = I.D. gripping	The base jaws move from the inside to the outside. The workpiece is gripped from the inside. With this setting, the gripper moves the base jaws together when released.
	FALSE = O.D. gripping	The base jaws move from the outside to the inside. The workpiece is gripped from the outside. With this setting, the gripper moves the base jaws apart when released.
Force	Value between 1 and 3	Adjustment of the gripping force in four steps: 1 = 100% 2 = 75% 3 = 50%
Target position	Value between 0.0 and max. stroke	The target position indicates up to where the gripper is to be opened or closed.

### 6.2.3 EGH Release – V100\_RELEASE

The gripper opens and releases the workpiece. The direction in which the gripper opens depends on the previous gripper command (V100\_GRIP).

### 6.2.4 EGH positioning movement – V100\_POS

This action allows the user to position the gripper fingers. The position can be absolute or relative.

#### CAUTION

#### Material damage to the workpiece possible!

This action is not suitable for gripping, as the movement always takes place at maximum speed.

- For positioning on an end stop, select the actions "Grip" and "Release".

#### NOTE

If the action "EGH positioning movement" is interrupted by a blockage, the drive switches off and an error message requiring acknowledgment is generated. The actuator remains de-energized until the next motion command.

#### Parameter settings V100\_POS

Gripping parameters	Value	Description
Target position	Value between 0.0 and max. stroke	The base jaws move to the position defined under "Target position".
isRelativeRun	FALSE = Positioning movement	The base jaws move to the position defined under "Target position".
	TRUE = Relative movement	The base jaws move to a position determined by the relative distance to the current position of the base jaws. The relative distance is also defined in the "Target position" field.
Wait	TRUE = Wait until position reached	If the command is selected, the entire movement waits before the next action in the robot program is executed.
	FALSE = Do not wait	If the command is not selected, the next component is executed during the position movement. This allows the robot arm to be moved during the position movement. No other gripping commands should follow in direct succession. Otherwise, the corresponding follow-up gripping command is executed before the movement is completed.

## 6.2.5 EGH light band - V100\_LIGHT

The signals of the light band are controlled by this component. The blue and red cable must be connected to the control box so that the light band can be controlled.

### Parameter settings V100\_LIGHT

Gripping parameters	Value	Description
deactivate	FALSE = LEDs on TRUE = LEDs off	Activates or deactivates the light band
setColor	[1] = green [2] = yellow [3] = red	Controls the colors on the light band

### Standard port assignment

- Blue cable = port 0 (P00)
- Red cable = port 1 (P01)

### Changing the port assignment


When inserting the component into the program sequence, two variables are automatically created. These define which cable is connected to which port.

- var\_blue\_wire\_EGH\_Lightband
- var\_red\_wire\_EGH\_Lightband

1. Open the "Variables" window in the main menu.

2. Change the ports.

- ⇒ Use values between 0 and 15 for both variables/port assignments (P00 – P015). The two variables may not have the same value.

Changing the port assignment of the light band

### 6.3 Using multiple grippers


Several grippers can be operated simultaneously on the robot via IO-Link. This requires adding a new device and changing the IO-Link ports.

The respective device number is displayed in the component name.

Component version	Function of the components	Gripper type	Device number
V100	RELEASE	EGH80	1

Tab.: Example name composition "V100\_RELEASE\_EGH801"

The components for the different devices are displayed in the program.


Display of two different devices in the program

#### Adding a new device

#### NOTE

The following applies in general: Only add as many devices as are actually connected to the robot.


1. Make sure that the desired component is already available in the program sequence (component for gripper 1).
2. Drag another component for gripper 2 (unit 2) into the sequence.
  - ⇒ A warning appears.
3. Select the "Add" button.

*Adding a new device*

⇒ If a second device already exists, it is displayed.

4. Select the desired device.
5. Select "Inherit old component" button to overwrite the component.


*Selecting an existing device*

⇒ A new device was added.

### Changing the IO-Link port

The communication parameters must be set. The procedure is identical for all components of the gripper except the light band component. This creates additional variables for setting the ports.


1. Open Modbus settings via  on the right side in the project window.
2. Select the component whose IO-Link port must be changed.
3. Select the "Import from IODD" button.

Modbus settings

**4.** Enter values:

- ⇒ IODD file: Select SCHUNK IODD file
- ⇒ Slave ID: 1
- ⇒ Enter start addresses for IO-Link ports according to the table

Entering values


Up to four grippers can be connected to the device. For an example see section ▶ 7.2 [ 66].

## 7 Example programs

### 7.1 Example: Initializing gripper

The following example shows how the gripper must be initialized at startup.

1. Acknowledge gripper with component "ACK".
  - ⇒ Status LEDs on the gripper change from red to green.
2. Perform reference run using the "REF" component.

Example - Initializing gripper

⇒ Gripper is ready for operation. Additional components can be used, ▶ 6.2 [58].


## 7.2 Example: Using multiple grippers

The following example shows how to initialize and change the IO-Link ports when using two grippers.

### Step 1: Initializing

■ The IODD file for the gripper is imported, ▶ 3.2 [44].

1. Initialize gripper 1. To do this, drag components "ACK" and "REF" into the program sequence.


#### Initializing gripper 1

2. Initialize gripper 2. To do this, drag components "ACK" and "REF" into the program sequence.


⇒ Warning message appears because the components "ACK" and "REF" for gripper 1 already exist.

3. Select the "Add" button.

#### Adding device 2


⇒ Both grippers are in the program sequence.

Gripper 1 and gripper 2 in the program sequence


## Step 2: Set IO-Link port "2"

1. Open Modbus settings via  on the right side in the project window.
2. First select component for gripper 2 and select "Edit" icon.

Gripper 2 - Modbus device settings


3. Select the "Import from IODD" button.

*Modbus settings*

**4. Enter values:**


- ⇒ IODD file: Select SCHUNK IODD file
- ⇒ Slave ID: 1
- ⇒ Start address in: Select IN address for port 2 (18)
- ⇒ Start address out: Select OUT address for port 2 (2065)



*Enter values*

**5. Set IO-Link port 2 for second component to analog.**

- ⇒ Two grippers can be controlled with the robot.

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