

Assembly and Operating Manual

SWM-B

Storage system

Translation of Original Operating
Manual

Imprint

Copyright:

This manual is protected by copyright. The author is SCHUNK SE & Co. KG.
All rights reserved.

Technical changes:

We reserve the right to make alterations for the purpose of technical improvement.

Document number: 1621461

Version: 03.00 | 31/01/2025 | en

Dear Customer,

Thank you for trusting our products and our family-owned company, the leading technology supplier of robots and production machines.

Our team is always available to answer any questions on this product and other solutions. Ask us questions and challenge us. We will find a solution!

Best regards,

Your SCHUNK team

Customer Management

Tel. +49-7133-103-2503

Fax +49-7133-103-2189

cmg@de.schunk.com



Please read the operating manual in full and keep it close to the product.

Table of Contents

1 General	4
1.1 About this manual.....	4
1.1.1 Applicable documents	4
1.1.2 Sizes.....	4
1.2 Warranty	5
1.3 Scope of delivery.....	5
1.4 Accessories	5
2 Basic safety notes	6
2.1 Appropriate use	6
2.2 Personal protective equipment	6
2.3 Constructional changes.....	6
2.4 Notes for assembly	6
2.5 Notes for operation	7
3 Technical data	8
4 Structure	9
5 Assembly and settings	10
5.1 Mounting the sensor (optional)	10
5.2 Mounting the storage module on the SCHUNK mounting plate.....	11
5.3 Mounting the storage module without the SCHUNK mounting plate.....	12
5.4 Mounting the adapter plate on the change system	13
5.5 Teaching routine	14
6 Translation of the original declaration of incorporation	16
7 Information on the RoHS Directive, REACH Regulation and Substances of Very High Concern (SVHC)	17

1 General

1.1 About this manual

This manual contains important information for a safe and appropriate use of the product.

This manual is an integral part of the product and must be kept accessible for the personnel at all times.

Before starting work, the personnel must have read and understood this operating manual. Prerequisite for safe working is the observance of all safety instructions in this manual.

In addition to these instructions, the documents listed under ► 1.1.1 [4] are applicable.

NOTE: The illustrations in this manual are intended to provide a basic understanding and may deviate from the actual version.

1.1.1 Applicable documents

- General terms of business *
- Catalog data sheet of the purchased product *

The documents labeled with an asterisk (*) can be downloaded from schunk.com/downloads.

1.1.2 Sizes

This operating manual applies to the following sizes:

- SWM-B 050
- SWM-B 085

1.2 Warranty

If the product is used as intended, the warranty is valid for 24 months from the ex-works delivery date under the following conditions:

- Observe the ambient conditions and operating conditions

Parts touching the workpiece and wear parts are not included in the warranty.

1.3 Scope of delivery

The scope of delivery includes

- Storage system SWM-B in the ordered size
- Assembly and Operating Manual
- Accessory pack

Content of the accessory pack:

- SWM-B 050: 2x Centering sleeve \emptyset 8 x 5.35
- SWM-B 085: 2x Centering sleeve \emptyset 12 x 6.65

1.4 Accessories

The following accessories are available for the product, which must be ordered separately:

- Adapter plates for quick-change adapter SWA
- Inductive proximity switch IN 5-S, Ident number: M8: 301469, M12: 301569
- Mounting plate for SWM-B, Ident number: 1547060

2 Basic safety notes

2.1 Appropriate use

- The SWM-B tool storage system is used for the safe storage of robot tools to enable an automated change of tools. For this purpose, the robot is equipped with a locking tool change system (e.g. SWS).

With the SWM-B, the tool is safely stored in a predefined orientation by its own weight and anti-rotation protection.

- The product may only be used within the scope of its technical data, ▶ 3 [8].
- The product is intended for industrial and industry-oriented use.
- The product is intended for installation in a machine/ automated system. The applicable guidelines for the machine/ automated system must be observed and complied with.
- Appropriate use of the product includes compliance with all instructions in this manual.

2.2 Personal protective equipment

- When working on and with the product, observe the occupational health and safety regulations and wear the required personal protective equipment.
- Wear protective gloves to guard against sharp edges and corners or rough surfaces.

2.3 Constructional changes

- Constructional changes may only be done with the permission of SCHUNK.

2.4 Notes for assembly

- Before starting assembly, secure the danger zone with suitable protective measures.

2.5 Notes for operation

- Observe safety distances.
- Never deactivate safety installations.
- If the energy supply is connected, do not move any parts by hand.
- Observe the valid country-specific safety and accident prevention regulations.
- Check the storage process after the robot crashes.

Possible electrostatic energy


Components or assembly groups may become electrostatically charged. When the electrostatic charge is touched, the discharge may trigger a shock reaction leading to injuries.

- The operator must ensure that all components and assembly groups are included in the local potential equalisation in accordance with the applicable regulations.
- While paying attention to the actual conditions of the working environment, the potential equalisation must be implemented by a specialist electrician according to the applicable regulations.
- The effectiveness of the potential equalisation must be verified by executing regular safety measurements.

3 Technical data


	SWM-B 050	SWM-B 085
Operating temperature [°C]		
Min.	5	5
Max.	60	60
Maximum permissible payload [kg]	16	100
Maximum permissible moment [Nm]		
M_x *	12	100
M_y	12	100

* Depending on the center of gravity, see diagrams below.



The tool center of gravity must be in the green area. In addition, the occurring moments must not exceed the permissible values.

4 Structure




- 1 Mounting plate (optional)
- 2 Storage module
- 3 Sensor for monitoring tool presence (optional)
- 4 Storage pins
- 5 Adapter plate A-SWA (in the illustration: A-SWA-007)

5 Assembly and settings

NOTE

The storage module may only be mounted in a horizontal position.


5.1 Mounting the sensor (optional)



The mounting screws for mounting the sensor are included in the scope of delivery of the sensor.

The cable outlet of the sensor is possible on both sides.

5.2 Mounting the storage module on the SCHUNK mounting plate



SWM-B 050 SWM-B 085

Mounting the mounting plate

Mounting screw (2x)	M8	M8
Mounting screw according to standard	ISO 4762	ISO 4762
Mounting screw strength class	8.8	8.8
Tightening torque [Nm]	20	20

Mounting the storage module on the mounting plate

Centering sleeve	∅ 8 x 5.35	∅ 12 x 6.65
Mounting screw (2x)	M4x20	M6x30
Mounting screw according to standard	DIN 6912 *	ISO 4762
Mounting screw strength class	8.8	12.9
Max. tightening torque [Nm]	3	15

* Alternatively, cylindrical screws according to ISO4762 can be used. These protrude beyond the housing, but this does not impair the function.

5.3 Mounting the storage module without the SCHUNK mounting plate

	SWM-B 050	SWM-B 085
Centering sleeve	∅ 8 x 5.35	∅ 12 x 5.65

Mounting the storage module from the front

Mounting screws	M4	M6
Mounting screw according to standard	DIN 6912	ISO 4762
Mounting screw strength class	8.8	12.9
Tightening torque [Nm]	3	15

Mounting the storage module from the rear


Mounting screws	M5	M8
Mounting screw strength class	8.8	8.8
Tightening torque [Nm]	5.5	20

The mounting point for the storage system must have a recess for the storage pin in the middle between the two mounting holes, as the storage pin is moved beyond the contour of the storage module during the storage process, see section "A" ▶ 5.5 [14].

	SWM-B 050	SWM-B 085
Depth [mm]	3	3
Diameter [mm]	10	15

Tab.: Recess dimensions

5.4 Mounting the adapter plate on the change system




	SWM-B 050	SWM-B 085
1	Side plate	
2	Intermediate plate	
3	Tightening torque of storage pin [Nm]	13
4	Mounting screws of the strength class	8.8

5.5 Teaching routine

Routine for removing the tool

1. Manually place the SWA quick-change adapter in the storage system.
2. Connect and lock the SWK quick-change master to the SWA quick-change adapter.



3. TCP programming.
The coordinates of the "TCP SWM-B" must be specified relative to the "TCP robot" (see tables).
4. Turn the tool with the robot by 10° to the horizontal axis around the "TCP SWM-B".
5. Move the tool linearly along the tilted X axis.

Ident number Intermediate plate	Size SWS	Z [mm]	X [mm]
1523807	-005	54.1	44.3
1523807	-011	46.7	44.3
1523816	-020	57.5	87.6
1523816	-021	61.0	87.6
1523816	-022	76.0	87.6
1523817	-040	77.3	100.6
1523819	-041	74.6	92.6
1523828	-060	62.3	97.6
1523830	-071	77.4	107.6

Ident number Intermediate plate	Size SWS	Z [mm]	X [mm]
1523831	-076	83.3	127.6

Tab.: TCP SWM-B coordinates for SWM-B variant with intermediate plate

Ident number Side plate	Size SWS	Z [mm]	X [mm]
1523812	-007	37.8	56.7
1523836	-029	51.0	96.6
1523836	-046	54.8	114.1

Tab.: TCP SWM-B coordinates for SWM-B variant with side plate

The values only apply to mounting the SWK quick-change master directly on the robot flange. Possible intermediate structures must also be taken into account.

Routine for storing the tool

To store the tool, carry out the described points in reverse order.

6 Translation of the original declaration of incorporation

in terms of the Directive 2006/42/EG, Annex II, Part 1 Section B.

Manufacturer/
Distributor SCHUNK SE & Co. KG
Spanntechnik | Greiftechnik | Automatisierungstechnik
Bahnhofstr. 106 – 134
D-74348 Lauffen/Neckar

We hereby declare that the partly completed machine described below

Product designation: Storage system / SWM-B /
ID number 1459336, 1459339

meets the following basic occupational health and safety of the Machinery Directive 2006/42/EC:

No. 1.1.1, No. 1.1.2, No. 1.1.3, No. 1.1.5, No. 1.3.2, No. 1.5.3, No. 1.5.4, No. 1.5.6, No. 1.5.8, No. 1.5.10, No. 1.5.11, No. 1.5.13

The partly completed machinery may not be put into operation until it has been confirmed that the machine into which the partly completed machinery is to be installed complies with the provisions of the Machinery Directive (2006/42/EC). The declaration shall be rendered invalid if modifications are made to the product.

Applied harmonized standards, especially:

EN ISO 12100:2010 Safety of machinery – General principles for design –
Risk assessment and risk reduction

The special technical documentation according to Annex VII, Part B, belonging to the partly completed machine, has been created.

Person authorized to compile the technical documentation:
Stefanie Walter, Address: see manufacturer's address

Signature: see original declaration

Lauffen/Neckar, January 2025

Dr.-Ing. Manuel Baumeister,
Head of Systems Engineering,
Technology & Innovation

7 Information on the RoHS Directive, REACH Regulation and Substances of Very High Concern (SVHC)

RoHS Directive

SCHUNK products are classified as "large-scale stationary installations" or as "large-scale stationary industrial tools" within the meaning of Directive 2011/65/EU and its extension 2015/863/EU "on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)", or fulfill their intended function only as part of one. Therefore products from SCHUNK do not fall within the scope of the directive at this time.

REACH Regulation


Products from SCHUNK fully comply with the regulations of Regulation (EC) No. 1907/2006 "concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH)" and its amendment 2022/477. SCHUNK attaches great importance to completely avoiding chemicals of concern to humans and the environment wherever possible.


Only in rare exceptional cases do SCHUNK products contain SVHC substances on the candidate list with a mass content above 0.1%. In accordance with Article. 33 (1) of Regulation (EC) No. 1907/2006, SCHUNK complies with its duty to "communicate information on substances in articles" and lists the components concerned and the substances used in an overview that can be viewed at [schunk.com/SVHC](https://www.schunk.com/SVHC).

Signature: see original declaration

Lauffen/Neckar, January 2025

Dr.-Ing. Manuel Baumeister,
Head of Systems Engineering,
Technology & Innovation







SCHUNK SE & Co. KG
Spanntechnik | Greiftechnik | Automatisierungstechnik

Bahnhofstr. 106 - 134
D-74348 Lauffen/Neckar
Tel. +49-7133-103-0
info@de.schunk.com
schunk.com

Folgen Sie uns | *Follow us*



Wir drucken nachhaltig | *We print sustainable*